

# Adapteur Power System II Console



*APS II Control Console  
AR-8300*

## APS II Control Console

This unique, versatile system is designed for multiple orthopaedic procedures. The APS II Control Console features two fully functional channels that allow simultaneous handpiece operation. Both channels can be operated with a single footswitch or individually, when a second footswitch is attached to the console. Additionally, either channel may be operated manually using the controls on the console face.

When simultaneously using both channels of the APS II console, each handpiece may be operated at the user's preferred mode and speed. Use of the gas pedal footswitch will override these settings, providing even greater flexibility in speed control.

### APS II Handpiece Speeds:

<b>Optional Handpiece Attachment</b>	<b>Speed</b>
Shaver Handpiece:	Forward/Reverse: Up to 8,000 rpm Oscillate: Up to 3,000 rpm
Drill Handpiece:	Forward/Reverse: Up to 1,400 rpm
Sagittal Saw Handpiece:	Forward/Reverse: 18,000 cpm

### APS II Console Specifications

Dimensions:	14.5" W x 5.2" H x 12.8" D (36.8 cm x 13.2 cm x 32.5 cm)
Power Requirements:	110 - 240 VAC, 6.3 A, 50/60 Hz, T 6.3 Amp
Power Consumption:	10 A / 250 V
Weight:	16.6 lbs (7.2 kg)